

TMEM91 Antibody



PACO30658

Product Information

Size:	Protein Background:
50ug	
Reactivity:	Gene ID:
Human	TMEM91
Source:	Uniprot
Rabbit	Q6ZNR0
Isotype:	Synonyms:
IgG	Transmembrane protein 91 (Dispanin subfamily C member 3) (DSPC3), TMEM91
Applications:	Immunogen:
ELISA, WB, IHC, IF	Recombinant Human Transmembrane protein 91 protein (1-97AA).
Recommended dilutions:	Storage:
ELISA:1:2000-1:10000, WB:1:1000-1:5000, IHC:1:20-1:200, IF:1:50-1:200	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western Blot

All lanes: TMEM91 antibody at 16 μ g/ml + A549 whole cell lysate


Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution


Predicted band size: 19, 16, 14, 15 kDa

Observed band size: 19, 34 kDa

.



Immunohistochemistry of paraffin-embedded human small intestine tissue using PACO30658 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human liver cancer using PACO30658 at dilution of 1:100.