

TMSB4Y Antibody



PACO30682

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Plays an important role in the organization of the cytoskeleton. Binds to and sequesters actin monomers (G actin) and therefore inhibits actin polymerization.

Gene ID:

TMSB4Y

Uniprot

O14604

Synonyms:

Thymosin beta-4, Y-chromosomal, TMSB4Y, TB4Y

Immunogen:

Recombinant Human Thymosin beta-4, Y-chromosomal protein (2-44AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A