

TNFAIP6 Antibody, HRP conjugated



PACO30695

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Possibly involved in cell-cell and cell-matrix interactions during inflammation and tumorigenesis.

Gene ID:

TNFAIP6

Uniprot

P98066

Synonyms:

Tumor necrosis factor-inducible gene 6 protein (Hyaluronate-binding protein) (TNF-stimulated gene 6 protein) (TSG-6) (Tumor necrosis factor alpha-induced protein 6) (TNF alpha-induced protein 6), TNFAIP6, TSG6

Immunogen:

Recombinant Human Tumor necrosis factor-inducible gene 6 protein (18-277AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A