

Tp53 Antibody



PACO30726

Product Information

Size:

50ug

Reactivity:

Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

Tp53

Uniprot

P02340

Synonyms:

Cellular tumor antigen p53 (Tumor suppressor p53), Tp53, P53 Trp53

Immunogen:

Recombinant Mouse Cellular tumor antigen p53 protein (1-387AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A