

PACO30730

Product Information

Size:

50ug

Reactivity:

Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

May function as an inhibitor of Wnt/beta-catenin signaling by indirectly interacting with LRP6 and blocking Wnt3a-dependent LRP6 internalization.

Gene ID:

Tpbg

Uniprot

Q5PQV5

Synonyms:

Trophoblast glycoprotein (5T4 oncofetal trophoblast glycoprotein) (5T4 oncotrophoblast glycoprotein) (Wnt-activated inhibitory factor 1) (WAI1), Tpbg, 5t4

Immunogen:

Recombinant Rat Trophoblast glycoprotein protein (32-361AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A