

TPPP2 Antibody, HRP conjugated



PACO30735

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

May bind tubulin but has no microtubule bundling activity.

Gene ID:

TPPP2

Uniprot

P59282

Synonyms:

Tubulin polymerization-promoting protein family member 2 (Protein p25-beta) (TPPP/p18), TPPP2, C14orf8

Immunogen:

Recombinant Human Tubulin polymerization-promoting protein family member 2 protein (1-170AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A