

TPPP3 Antibody



PACO30738

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Binds tubulin and has microtubule bundling activity. May play a role in cell proliferation and mitosis.

Gene ID:

TPPP3

Uniprot

Q9BW30

Synonyms:

Tubulin polymerization-promoting protein family member 3 (TPPP/p20), TPPP3


Immunogen:

Recombinant Human Tubulin polymerization-promoting protein family member 3 protein (2-176AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human testis tissue using PACO30738 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human kidney tissue using PACO30738 at dilution of 1:100.