

# TRAF5 Antibody, FITC conjugated



PACO30776

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Adapter protein and signal transducer that links members of the tumor necrosis factor receptor family to different signaling pathways by association with the receptor cytoplasmic domain and kinases. Mediates activation of NF-kappa-B and probably JNK. Seems to be involved in apoptosis. Plays a role in mediating activation of NF-kappa-B by EIF2AK2/PKR.

**Gene ID:**

TRAF5

**Uniprot**

O00463

**Synonyms:**

TNF receptor-associated factor 5 (RING finger protein 84), TRAF5, RNF84

**Immunogen:**

Recombinant Human TNF receptor-associated factor 5 protein (1-557AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A