

TRAF5 Antibody, Biotin conjugated



PACO30777

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Adapter protein and signal transducer that links members of the tumor necrosis factor receptor family to different signaling pathways by association with the receptor cytoplasmic domain and kinases. Mediates activation of NF-kappa-B and probably JNK. Seems to be involved in apoptosis. Plays a role in mediating activation of NF-kappa-B by EIF2AK2/PKR.

Gene ID:

TRAF5

Uniprot

O00463

Synonyms:

TNF receptor-associated factor 5 (RING finger protein 84), TRAF5, RNF84

Immunogen:

Recombinant Human TNF receptor-associated factor 5 protein (1-557AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A