

RPL17 Antibody



PACO30830

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:5000

Protein Background:

RPL17 is a ribosomal protein that is a component of the large 60S subunit. The protein belongs to the L22P family of ribosomal proteins.

Gene ID:

RPL17

Uniprot

P18621

Synonyms:

60S ribosomal protein L17 (60S ribosomal protein L23) (Large ribosomal subunit protein uL22) (PD-1), RPL17

Immunogen:

Recombinant Human 60S ribosomal protein L17 protein (2-184AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western blot

All lanes: RPL17 antibody at 2 μ g/ml + EC109 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/15000 dilution

Predicted band size: 22, 18, 27 kDa

Observed band size: 27, 80 kDa