

# RPL17 Antibody, HRP conjugated



PACO30835

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

RPL17 is a ribosomal protein that is a component of the large 60S subunit. The protein belongs to the L22P family of ribosomal proteins.

**Gene ID:**

RPL17

**Uniprot**

P18621

**Synonyms:**

60S ribosomal protein L17 (60S ribosomal protein L23) (Large ribosomal subunit protein uL22) (PD-1), RPL17

**Immunogen:**

Recombinant Human 60S ribosomal protein L17 protein (2-184AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A