

TSLP Antibody, Biotin conjugated



PACO30841

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Isoform 1: Cytokine that induces the release of T-cell-attracting chemokines from monocytes and, in particular, enhances the maturation of CD11c(+) dendritic cells. Can induce allergic inflammation by directly activating mast cells.

Gene ID:

TSLP

Uniprot

Q969D9

Synonyms:

Thymic stromal lymphopoietin, TSLP

Immunogen:

Recombinant Human Thymic stromal lymphopoietin protein (29-159AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A