

TSPAN1 Antibody



PACO30842

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200,
IF:1:50-1:200

Protein Background:

cytoplasm, extracellular exosome, integral component of plasma membrane, membrane, perinuclear region of cytoplasm, plasma membrane, vesicle, cell migration, cell proliferation, cell surface receptor signaling pathway.

Gene ID:

TSPAN1

Uniprot

O60635

Synonyms:

Tetraspanin-1 (Tspan-1) (Tetraspan NET-1) (Tetraspanin TM4-C), TSPAN1


Immunogen:

Recombinant Human Tetraspanin-1 protein (110-211AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human colon cancer using PACO30842 at dilution of 1:100.



Immunofluorescent analysis of HeLa cells using PACO30842 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).