

# Ttr Antibody, Biotin conjugated



PACO30884

---

## Product Information

**Size:**

50ug

**Reactivity:**

Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Thyroid hormone-binding protein. Probably transports thyroxine from the bloodstream to the brain.

**Gene ID:**

Ttr

**Uniprot**

P07309

**Synonyms:**

Transthyretin (Prealbumin), Ttr

**Immunogen:**

Recombinant Mouse Transthyretin protein (23-147AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A