
Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Activates ubiquitin by first adenylating with ATP its C-terminal glycine residue and thereafter linking this residue to the side chain of a cysteine residue in E1, yielding a ubiquitin-E1 thioester and free AMP. Catalyzes the ISGylation of influenza A virus NS1 protein.

Gene ID:

UBA7

Uniprot

P41226

Synonyms:

Ubiquitin-like modifier-activating enzyme 7 (Ubiquitin-activating enzyme 7) (D8)
(Ubiquitin-activating enzyme E1 homolog), UBA7, UBE1L UBE2

Immunogen:

Recombinant Human Ubiquitin-like modifier-activating enzyme 7 protein (1-265AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A