

# UCP2 Antibody, Biotin conjugated



PACO30983

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

UCP are mitochondrial transporter proteins that create proton leaks across the inner mitochondrial membrane, thus uncoupling oxidative phosphorylation from ATP synthesis. As a result, energy is dissipated in the form of heat.

**Gene ID:**

UCP2

**Uniprot**

P55851

**Synonyms:**

Mitochondrial uncoupling protein 2 (UCP 2) (Solute carrier family 25 member 8) (UCPH), UCP2, SLC25A8

**Immunogen:**

Recombinant Human Mitochondrial uncoupling protein 2 protein (137-180AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A