

UHRF1BP1L Antibody



PACO30996

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

UHRF1BP1L

Uniprot

A0JNW5

Synonyms:

UHRF1-binding protein 1-like, UHRF1BP1L, KIAA0701

Immunogen:

Recombinant Human UHRF1-binding protein 1-like protein (251-522AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A