

PACO31032

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

May participate in trafficking events that are associated with myogenesis, such as myoblast fusion and/or GLUT4 trafficking.

Gene ID:

VAMP5

Uniprot

O95183

Synonyms:

Vesicle-associated membrane protein 5 (VAMP-5) (Myobrevin), VAMP5


Immunogen:

Recombinant Human Vesicle-associated membrane protein 5 protein (1-72AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human lung cancer using PACO31032 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human colon cancer using PACO31032 at dilution of 1:100.