

ZNF155 Antibody



PACO31160

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

May be involved in transcriptional regulation.

Gene ID:

ZNF155

Uniprot

Q12901

Synonyms:

Zinc finger protein 155, ZNF155


Immunogen:

Recombinant Human Zinc finger protein 155 protein (1-538AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human liver cancer using PACO31160 at dilution of 1:100.