

# ZNF165 Antibody



PACO31164

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200

**Protein Background:**

May be involved in transcriptional regulation.

**Gene ID:**

ZNF165

**Uniprot**

P49910

**Synonyms:**

Zinc finger protein 165 (Cancer/testis antigen 53) (CT53) (LD65) (Zinc finger and SCAN domain-containing protein 7), ZNF165, ZPF165 ZSCAN7

**Immunogen:**


Recombinant Human Zinc finger protein 165 protein (1-485AA).

**Storage:**


Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---



Immunohistochemistry of paraffin-embedded human testis tissue using PACO31164 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human cervical cancer using PACO31164 at dilution of 1:100.