

ZNF19 Antibody



PACO31188

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:2000,
IHC:1:20-1:200

Protein Background:

May be involved in transcriptional regulation.

Gene ID:

ZNF19

Uniprot

P17023

Synonyms:

Zinc finger protein 19 (Zinc finger protein KOX12), ZNF19, KOX12


Immunogen:

Recombinant Human Zinc finger protein 19 protein (1-200AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western blot

All lanes: ZNF19 antibody at 4 μ g/ml

Lane 1: Mouse liver tissue

Lane 2: MDA-MB-231 whole cell lysate

Lane 3: HepG2 whole cell lysate


Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution


Predicted band size: 53, 52 kDa

Observed band size: 53 kDa

.



Immunohistochemistry of paraffin-embedded human small intestine tissue using PACO31188 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human breast cancer tissue using PACO31188 at dilution of 1:100.