

# ZNF212 Antibody



PACO31200

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, WB:1:500-1:2000,  
IHC:1:20-1:200

**Protein Background:**

May be involved in transcriptional regulation.

**Gene ID:**

ZNF212

**Uniprot**

Q9UDV6

**Synonyms:**

Zinc finger protein 212 (Zinc finger protein C2H2-150), ZNF212, ZNFC150

**Immunogen:**


Recombinant Human Zinc finger protein 212 protein (1-495AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---



### Western blot


All lanes: ZNF212 antibody at 8 $\mu$ g/ml + Mouse liver tissue

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 56 kDa

Observed band size: 56 kDa



Immunohistochemistry of paraffin-embedded human lung tissue using PACO31200 at dilution of 1:100.