

ZNF223 Antibody



PACO31212

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

May be involved in transcriptional regulation.

Gene ID:

ZNF223

Uniprot

Q9UK11

Synonyms:

Zinc finger protein 223, ZNF223

Immunogen:

Recombinant Human Zinc finger protein 223 protein (201-482AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A