

ZNF287 Antibody



PACO31252

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

May be involved in transcriptional regulation.

Gene ID:

ZNF287

Uniprot

Q9HBT7

Synonyms:

Zinc finger protein 287 (Zinc finger protein with KRAB and SCAN domains 13), ZNF287, ZKSCAN13


Immunogen:

Recombinant Human Zinc finger protein 287 protein (130-364AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human pancreatic cancer using PACO31252 at dilution of 1:100.