

ZNF32 Antibody



PACO31260

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

May be involved in transcriptional regulation.

Gene ID:

ZNF32

Uniprot

P17041

Synonyms:

Zinc finger protein 32 (C2H2-546) (Zinc finger protein KOX30), ZNF32, KOX30

Immunogen:

Recombinant Human Zinc finger protein 32 protein (1-273AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A