

ZNF385C Antibody



PACO31280

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

p53 binding, RNA binding, intrinsic apoptotic signaling pathway by p53 class mediator.

Gene ID:

ZNF385C

Uniprot

Q66K41

Synonyms:

Zinc finger protein 385C, ZNF385C

Immunogen:

Recombinant Human Zinc finger protein 385C protein (1-174AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A