

ZNF397 Antibody



PACO31288

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Isoform 3 acts as a DNA-dependent transcriptional repressor.

Gene ID:

ZNF397

Uniprot

Q8NF99

Synonyms:

Zinc finger protein 397 (Zinc finger and SCAN domain-containing protein 15) (Zinc finger protein 47), ZNF397, ZNF47 ZSCAN15


Immunogen:

Recombinant Human Zinc finger protein 397 protein (1-275AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human testis tissue using PACO31288 at dilution of 1:100.