

# ZNF433 Antibody



PACO31308

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

May be involved in transcriptional regulation.

**Gene ID:**

ZNF433

**Uniprot**

Q8N7K0

**Synonyms:**

Zinc finger protein 433, ZNF433

**Immunogen:**

Recombinant Human Zinc finger protein 433 protein (1-300AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A