

# ZSCAN32 Antibody



PACO31312

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

ELISA:1:2000-1:10000, WB:1:500-1:2000

**Protein Background:**

May be involved in transcriptional regulation.

**Gene ID:**

ZSCAN32

**Uniprot**

Q9NX65

**Synonyms:**

Zinc finger and SCAN domain-containing protein 32 (Human cervical cancer suppressor gene 5 protein) (HCCS-5) (Zinc finger protein 434), ZSCAN32, ZNF434

**Immunogen:**

Recombinant Human Zinc finger and SCAN domain-containing protein 32 protein (6-248AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---



### Western blot

All lanes: ZSCAN32 antibody at 1 $\mu$ g/ml + MCF-7 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 79, 55, 47 kDa

Observed band size: 55 kDa

.