

ZNF490 Antibody



PACO31340

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

May be involved in transcriptional regulation.

Gene ID:

ZNF490

Uniprot

Q9ULM2

Synonyms:

Zinc finger protein 490, ZNF490, KIAA1198


Immunogen:

Recombinant Human Zinc finger protein 490 protein (101-400AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human testis tissue using PACO31340 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human colon cancer using PACO31340 at dilution of 1:100.