

ZNF513 Antibody



PACO31348

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Transcriptional regulator that plays a role in retinal development and maintenance.

Gene ID:

ZNF513

Uniprot

Q8N8E2

Synonyms:

Zinc finger protein 513, ZNF513


Immunogen:

Recombinant Human Zinc finger protein 513 protein (1-300AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human skin tissue using PACO31348 at dilution of 1:100.