

ZNF530 Antibody



PACO31356

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

May be involved in transcriptional regulation.

Gene ID:

ZNF530

Uniprot

Q6P9A1

Synonyms:

Zinc finger protein 530, ZNF530, KIAA1508

Immunogen:

Recombinant Human Zinc finger protein 530 protein (301-599AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A