

ZNF563 Antibody



PACO31384

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:1000-1:5000,
IHC:1:20-1:200

Protein Background:

May be involved in transcriptional regulation.

Gene ID:

ZNF563

Uniprot

Q8TA94

Synonyms:

Zinc finger protein 563, ZNF563


Immunogen:

Recombinant Human Zinc finger protein 563 protein (214-476AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western blot

All lanes: ZNF563 antibody at 7 μ g/ml

Lane 1: HeLa whole cell lysate

Lane 2: A549 whole cell lysate

Lane 3: MCF-7 whole cell lysate

Lane 4: K562 whole cell lysate


Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 56, 27 kDa

Observed band size: 56 kDa

.



Immunohistochemistry of paraffin-embedded human kidney tissue using PACO31384 at dilution of 1:100.