

# ZNF606 Antibody



PACO31412

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

May act as a transcriptional repressor.

**Gene ID:**

ZNF606

**Uniprot**

Q8WXB4

**Synonyms:**

Zinc finger protein 606 (Zinc finger protein 328), ZNF606, KIAA1852 ZNF328

**Immunogen:**

Recombinant Human Zinc finger protein 606 protein (501-792AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A