

# ZNF622 Antibody



PACO31428

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200

**Protein Background:**

May behave as an activator of the bound transcription factor, MYBL2, and be involved in embryonic development.

**Gene ID:**

ZNF622

**Uniprot**

Q969S3

**Synonyms:**

Zinc finger protein 622 (Zinc finger-like protein 9), ZNF622, ZPR9

**Immunogen:**


Recombinant Human Zinc finger protein 622 protein (2-280AA).

**Storage:**


Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---



Immunohistochemistry of paraffin-embedded human liver cancer using PACO31428 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human testis tissue using PACO31428 at dilution of 1:100.