

ZNF626 Antibody



PACO31440

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:2000,
IHC:1:20-1:200

Protein Background:

May be involved in transcriptional regulation.

Gene ID:

ZNF626

Uniprot

Q68DY1

Synonyms:

Zinc finger protein 626, ZNF626


Immunogen:

Recombinant Human Zinc finger protein 626 protein (1-92AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western blot

All lanes: ZNF626 antibody at 10 μ g/ml

Lane 1: HepG2 whole cell lysate

Lane 2: Raji whole cell lysate


Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 61, 53, 11 kDa

Observed band size: 61 kDa

.



Immunohistochemistry of paraffin-embedded human breast cancer using PACO31440 at dilution of 1:100.