

ZNF689 Antibody



PACO31472

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:1000-1:5000

Protein Background:

May be involved in transcriptional regulation.

Gene ID:

ZNF689

Uniprot

Q96CS4

Synonyms:

Zinc finger protein 689, ZNF689


Immunogen:

Recombinant Human Zinc finger protein 689 protein (1-300AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western blot

All lanes: ZNF689 antibody at 5 μ g/ml

Lane 1: Hela whole cell lysate

Lane 2: HepG2 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 57 kDa

Observed band size: 57 kDa

.