

# ZNF707 Antibody



PACO31480

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, WB:1:500-1:5000,  
IHC:1:20-1:200

**Protein Background:**

May be involved in transcriptional regulation.

**Gene ID:**

ZNF707

**Uniprot**

Q96C28

**Synonyms:**

Zinc finger protein 707, ZNF707

**Immunogen:**


Recombinant Human Zinc finger protein 707 protein (51-369AA).

**Storage:**


Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---



Immunohistochemistry of paraffin-embedded human kidney tissue using PACO31480 at dilution of 1:100.



### Western Blot

Positive WB detected in: Mouse heart tissue, Mouse stomach tissue, Mouse skeletal muscle tissue


All lanes: ZNF707 antibody at 3.4µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 44 kDa

Observed band size: 44 kDa



Immunohistochemistry of paraffin-embedded human gastric cancer tissue using PACO31480 at dilution of 1:100.