

ZNF75A Antibody



PACO31484

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:1000-1:5000

Protein Background:

May be involved in transcriptional regulation.

Gene ID:

ZNF75A

Uniprot

Q96N20

Synonyms:

Zinc finger protein 75A, ZNF75A


Immunogen:

Recombinant Human Zinc finger protein 75A protein (1-296AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western blot

All lanes: ZNF75A antibody at 10 μ g/ml

Lane 1: Mouse brain tissue

Lane 2: Mouse liver tissue

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 35 kDa

Observed band size: 35, 65 kDa