

ZNF829 Antibody



PACO31516

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:2000

Protein Background:

May be involved in transcriptional regulation.

Gene ID:

ZNF829

Uniprot

Q3KNS6

Synonyms:

Zinc finger protein 829, ZNF829


Immunogen:

Recombinant Human Zinc finger protein 829 protein (133-432AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western Blot

Positive WB detected in: HeLa whole cell lysate

All lanes: ZNF829 antibody at 3 μ g/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 51, 29, 59 kDa

Observed band size: 51 kDa