

ERH Antibody



PACO31544

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

May have a role in the cell cycle.

Gene ID:

ERH

Uniprot

P84090

Synonyms:

Enhancer of rudimentary homolog, ERH


Immunogen:

Recombinant Human Enhancer of rudimentary homolog protein (2-104AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human lung tissue using PACO31544 at dilution of 1:100.