

C1D Antibody, Biotin conjugated



PACO31691

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Plays a role in the recruitment of the RNA exosome complex to pre-rRNA to mediate the 3-5end processing of the 5.8S rRNA; this function may include MPHOSPH6. Can activate PRKDC not only in the presence of linear DNA but also in the presence of supercoiled DNA. Can induce apoptosis in a p53/TP53 dependent manner. May regulate the TRAX/TSN complex formation. Potentiates transcriptional repression by NR1D1 and THRB.

Gene ID:

C1D

Uniprot

Q13901

Synonyms:

Nuclear nucleic acid, binding protein C1D (hC1D), C1D

Immunogen:

Recombinant Human Nuclear nucleic acid, binding protein C1D protein (1-135AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A