

# TADA3 Antibody



PACO31700

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Functions as a component of the PCAF complex. The PCAF complex is capable of efficiently acetylating histones in a nucleosomal context. The PCAF complex could be considered as the human version of the yeast SAGA complex. Also known as a coactivator for p53/TP53-dependent transcriptional activation. Component of the ATAC complex, a complex with histone acetyltransferase activity on histones H3 and H4.

**Gene ID:**

TADA3

**Uniprot**

O75528

**Synonyms:**

Transcriptional adapter 3 (ADA3 homolog) (hADA3) (STAF54) (Transcriptional adapter 3-like) (ADA3-like protein), TADA3, ADA3 TADA3L

**Immunogen:**

Recombinant Human Transcriptional adapter 3 protein (1-238AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A