

ZNF592 Antibody



PACO31724

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

May be involved in transcriptional regulation.

Gene ID:

ZNF592

Uniprot

Q92610

Synonyms:

Zinc finger protein 592, ZNF592, KIAA0211

Immunogen:

Recombinant Human Zinc finger protein 592 protein (1-242AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A