

# SNX22 Antibody



PACO31800

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200

**Protein Background:**

May be involved in several stages of intracellular trafficking. Interacts with membranes containing phosphatidylinositol 3-phosphate (PtdIns(3P)).

**Gene ID:**

SNX22

**Uniprot**

Q96L94

**Synonyms:**

Sorting nexin-22, SNX22

**Immunogen:**


Recombinant Human Sorting nexin-22 protein (1-120AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---



Immunohistochemistry of paraffin-embedded human tonsil tissue using PACO31800 at dilution of 1:100.