

EFCAB14 Antibody



PACO31848

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:2000,
IHC:1:20-1:200

Protein Background:**Gene ID:**

EFCAB14

Uniprot

O75071

Synonyms:

EF-hand calcium-binding domain-containing protein 14, EFCAB14, KIAA0494


Immunogen:

Recombinant Human EF-hand calcium-binding domain-containing protein 14 protein (1-483AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western blot

All lanes: EFCAB14 antibody at 4 μ g/ml

Lane 1: Colo320 whole cell lysate

Lane 2: Hela whole cell lysate


Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 56 kDa

Observed band size: 56 kDa

.



Immunohistochemistry of paraffin-embedded human thyroid tissue using PACO31848 at dilution of 1:100.