

MRPL9 Antibody



PACO31876

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:5000

Protein Background:**Gene ID:**

MRPL9

Uniprot

Q9BYD2

Synonyms:

39S ribosomal protein L9, mitochondrial (L9mt) (MRP-L9) (Mitochondrial large ribosomal subunit protein bL9m), MRPL9


Immunogen:

Recombinant Human 39S ribosomal protein L9, mitochondrial protein (60-267AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western blot

All lanes: MRPL9 antibody at 2 μ g/ml

Lane 1: EC109 whole cell lysate

Lane 2: 293T whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/15000 dilution

Predicted band size: 31 kDa

Observed band size: 31 kDa

.