

IFNB1 Antibody



PACO32192

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Has antiviral, antibacterial and anticancer activities.

Gene ID:

IFNB1

Uniprot

P01574

Synonyms:

Interferon beta (IFN-beta) (Fibroblast interferon), IFNB1, IFB IFNB


Immunogen:

Recombinant Human Interferon β protein (22-187AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human liver cancer using PACO32192 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human colon cancer using PACO32192 at dilution of 1:100.