

MPZ Antibody, HRP conjugated



PACO32477

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

MPZ

Uniprot

P25189

Synonyms:

Myelin protein P0 (Myelin peripheral protein) (MPP) (Myelin protein zero), MPZ

Immunogen:

Recombinant Human Myelin protein P0 protein (30-156AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A