

PMP2 Antibody, HRP conjugated



PACO32481

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

PMP2 is a small basic protein, belonging to the fatty-acid, binding protein (FABP) family and is one of the major proteins of peripheral myelin. It appears to be involved in the transport of fatty acid, or the metabolism of myelin lipids, and has been shown to have lipid-binding activity. PMP2 protein may have an important role in the organization of compact myelin.

Gene ID:

PMP2

Uniprot

P02689

Synonyms:

Myelin P2 protein (Peripheral myelin protein 2), PMP2

Immunogen:

Recombinant Human Myelin P2 protein (1-132AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A